



The Midmark 411 is the ultimate universal power table because it offers everything physicians need for their medical environment. With its ergonomic design, the Midmark 411 not only looks great... but it comes standard with both a programmable foot and hand control for added convenience. The table can be programmed to four positions, and will move quickly and precisely where the physician wants it... every time.

Features

- The Midmark 411 Power Chair is a refurbished 4- way power chair designed for multi-practice situations. This chair operates with a side panel control to control the back, tilt, foot and height.
- The 411 is fully programmable with 8 program positions. This feature saves time during patient exams to get the table into and out of common positions.
- Other useful features include procto knee rest, paper storage, slide-out treatment pan, adjustable headrest and adjustable stirrups, removable upholstery sections for effective and easy cleaning, optional foot control.
- Perfect for plastic surgery, OB-Gyn, proctology, dermatology and general practice.



Specifications

Dimensions

Upholstered top: 27.0 in. wide x 69.25 in long (68.6 cm x 175.9 cm)
With Headrest Extended: 82.75 in (210.2 cm)

Specifications

Patient Weight (Maximum): 325 lbs (147.4 kgs) maximum
Weight of Table: 515 lbs (233.6 kgs)
Back Section Range: 0° (horizontal) to +85° ± 5°
Table Top Height Range: 22.5 ± 0.5 in. to 40.0 ± 1.0 in. (55.9 ± 1.3 cm to 101.6 ± 2.5 cm)
Foot Section Range: 0° (horizontal) to -90° ± 5°
Table Top Tilt Range: 0° (horizontal) to +30° ± 5° (foot up)

Electrical Requirments

Electrical: 230 VAC nominal (220 VAC minimum - 240 VAC maximum), 50/60 HZ, 8 amps (Max)
Paper Rolls: Can accept two standard paper rolls of 18 in. x 3.0 in. (45.7 cm x 7.6 cm)
Duty Cycle: 15 seconds on / 5 minutes off (motor run time)
Power Cord: Extends 70 in. (Minimum) from table. 12 AWG / 3 conductor, SJT grey jacketed junior hard service with hospital grade grounding type plug.