

## Midmark Ritter M11 Autoclaves



### Features:

Designed for fast, efficient sterilization of heat- and moisture-stable materials, prevacuum sterilizers are equipped with prevacuum, gravity, flash, express, leak test, and daily air removal test cycles.

### Specifications:

**Dimensions:** 24 x 36 x 48 inches (610 x 914 x 1219 MM)

**Pressure Chamber, Recessed Installation, Radial Arm Door, Printcon Controls, Integrated Thermal Array Recorder**



Typical manufacturer's specifications and images