

# Nellcor NBP-75 Monitor



Capnography For mechanically ventilated and spontaneously breathing patients. Full-featured capnography previously available only in much larger monitors. For use on intubated and nonintubated patients. Low sample flow rate of 50 ml/min allows monitoring on the widest range of patients, from neonates to adults. Patented Microstream™ technology and versatile range of FilterLine accessories enable CO<sub>2</sub> monitoring even in challenging high-humidity environments. Innovative optical bench enhances stable, accurate measurements from a small sample. Capnographic waveforms and trends. Pulse Oximetry NPB-75 Trusted Nellcor pulse oximetry technology. Compatible with complete line of Nellcor adhesive and reusable sensors. Graphical trends and plethysmographic waveform.

