



Features of the Nellcor NPB 70

- Full-featured capnography previously available only in much larger monitors.
- Innovative optical bench enhances stable, accurate measurements from a small sample.
- Low sample flow rate of 50 mL/min allows monitoring on wide patient range without compromising response time.
- Capnographic waveforms and trends.
- Data output/printing.



Nellcor NPB 70 Capnography Monitor

Indications for Usage

- For use on neonatal to adult patients
- Operates in high humidity environments
- Used in most any clinical settings
- For use with intubated or nonintubated patients

Physical Characteristics:

Size:	8.11 H x 3.46 W x 2.06 D inches (20.6 x 8.8 x 5.2 cm).
Weight (with batteries):	1.66 lb (750 g).



Soma Technology Inc.
166 Highland Park Drive
Bloomfield, CT 06002 • USA
1.800.438.7662
www.SomaTechnology.com

Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.