



## Features

- Easy to operate and maintain
- Single or dual endoscope reprocessing capability
- Minimizes liquid chemical germicide fume exposure
- Semi-automated alcohol rinsing and air drying of endoscope channels
- Heated liquid chemical germicide reservoir
- Fast disinfection cycle turnaround time with Rapicide™ by MediVators
- Compatible with Olympus® ultrasound endoscopes
- Disinfectant immersion times for Rapicide™, glutaraldehyde, and o-phthalaldehyde
- Two-rinse and three-rinse protocols available

