

Philips IntelliVue MP50




Features

- 30.5cm (12") color touchscreen monitor display with 4 or 6 waveforms it's bright and easy to read.
- Comfortable handle and rugged housing for easy portability.
- IntelliVue comes with 10 preset screen configurations.
- Customized viewing options let you view and analyze data in graphical or numerical formats, juxtapose real-time measurements and trended data, and organize every onscreen element – from waveforms to data labels – as desired.
- Touchscreen operation makes many functions directly accessible through simple commands.
- Access information from a variety of systems on the hospital network – pharmacy, laboratory, ECG management systems, radiology PACS, host based hospital systems and others – and deliver it to the patient monitor through IntelliVue's portal technology.
- No separate hard drive and no fan. IntelliVue is space-saving, quiet, stable, and starts quickly.
- Navigation Point operation designed for easy information input and onscreen navigation.
- Options to add Anesthesia, Intensive Care Package or Neonatal Monitoring



Note: The technical data given in this publication is for general information and is subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.



Soma Technology Inc. • 166 Highland Park Drive • Bloomfield, CT 06002 • USA
Phone: 1.800.GET.SOMA • www.SomaTechnology.com • Email: soma@somatechnology.com



Revised 5/15/2013