



Built to work for all types of general surgery, the Amsco 3080 operating table combines easy patient positioning, good C-arm access, and convenient usage into one reliable surgical table. The versatility and simplicity of the Amsco 3080 table provide excellent imaging capabilities during even the most demanding procedures. With override protection, the Amsco 3080 O.R. table can still position patients or return the table to the neutral position in the event of electrical malfunction or complete power failure.

Features

- Surgical Table with Radio Top
- One of the most durable, flexible and reliable surgical tables available
- Brings more efficiency to the operating room
- Fluoroscopic capabilities
- Simple patient positioning
- Returns to level with one simple command



Specifications

Dimensions

Length: 75.94".
Width: 20".
Width with Rail: 21.63"
Weight limit: 500 lbs

Table Height:

Raised: 44".
Lowered: 27".
Trendelenburg: 25°.
Reverse Trendelenburg: 25°.

Lateral Tilt

± 18°

Back Section

Down: 25°.
Up: 55°

Leg Section

Down: 105°.
Up: 80°.
Reflex: 20°.
Flex: 30°.

Kidney Bridge

4.5" elevation.

