



The **Steris Amsco Century V116** is a full-size steam sterilizer system. The Century V116 offers a fast and efficient sterilization process. The sterilizer comes with a touchscreen display. This makes it easy to select any one of the pre-programmed cycles, including; pre-vac, gravity, flash, express, leak test, and a daily air removal testing cycle. The Chamber is 16 inches x 16 inches x 20 inches and is opened through a power sliding door operated by a foot pedal. The Operation temperature of the sterilizer ranges from 250° to 270°. The Steris Amsco Century Series of Sterilizers also includes the Steris Amsco Century V120.

Features

- Cabinet Mounted with a Stainless Steel Enclosure
- Foot Pedal Actuated Door
- One (1) Rack
- Optional built-in steam generator
- Touch Screen Display
- Two (2) Shelves
- Vertical Door: Automatic



Specifications

Chamber Dimensions	16" x 16" x 26"
Overall Exterior Dimensions (including door lock-wheel)	74.5"H x 30"W x 35.75"D (with door closed)
Operating Temperature	250° to 270°

