



The Stryker LD304 Birthing Bed provides comfort and security for both the mother and caregivers. This birthing bed is simple to operate, mobile and it's designed to provide safety while bringing comfort to mothers in labor. Convenient controls provide the Stryker birthing bed to get into trend positions with a single button, these options making it easy to customize the positions for minimum effort.

The Stryker foot section promotes proper body mechanics to help prevent caregiver back strain or injury. The Stryker Birthing Bed features the patented Skoocher, which moves the patient or mother closer to the perineal edge of the bed with the touch of a button. This is a very secure and comfortable bed for mothers and caregivers alike, and its simple operation, mobility, and ergonomic features are all you're looking for in a birthing bed.

Features

- Lock-Rite slide-off foot section
- Lightweight, self-standing foot section
- Glideaway foot supports with six-way adjustability
- Single-button Trend.
- Three-position side rails
- In-rail motion controls with lockout
- 17.5-inch (44 cm) low bed height
- 5-inch (13 cm) Ultra Comfort head and a seat support surface
- 3-inch (8 cm) foot mattress
- 6-inch (15 cm) Omni Surface casters with wheel covers
- Locking-caster steering
- Four-wheel, steel-ring brake system with centrally located activation pedal
- 33-inch (84 cm) mattress width with four motor functions (head, bed, foot, Trend.)
- Automatic pelvic tilt
- CPR release
- Manual backup



Specifications

Overall length	93" (236 cm)
Overall Width	Siderails: Up 41" (104 cm) Siderails: Down 37" (94 cm)
Lift Capacity	Bed: 500 lbs (228 cm) Head Section: 250lbs (114 cm) Foot Section: 300lbs (136 cm)
Height Range (to litter top)	High: 35.5" (90 cm) Low: 17.5" (44 cm)
Litter Positioning	Backrest: 0° – 70° Pelvic Tilt: 4° Maximum Trend
Patient Surface	33" x 81" (84cm x 206 cm)
Mattress Thickness	Head/Seat: 5" (13 cm) Foot: 3" (8cm)
Seat Depth	Standard: 12.5" (32 cm) With Skoocher: 2.5" - 8.5" (32 - 22cm)
Caster Diameter	Standard: 6" (15 cm) Optional: 18" (20cm)