



Large Volume Steam Autoclave by Tuttnauer. Tuttnauer 3870EA Autoclave features a 15" Chamber. The 3870EA Autoclave is the largest of Tuttnauer's offering. This Autoclave features a closed door active drying system. Perfect for multi-doctor practices, clinics or surgery centers. 15"x30" chamber accommodates longer instruments or a large number of small instruments.

## Features

- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Overheat cutoff prevents burnout of heating elements in case of insufficient water in the chamber.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.
- Closed Door Dry Cycle
- Integrated printer for cycle documentation (optional)
- Chamber Dimensions: 15.1" x 29.9"
- Overall Dimensions: 34.5" x 26" x 20.6"
- Cycle Time - 55 min (cold), 40 min (hot)



# Specifications

## Technical Specifications

Chamber Dimensions x L - inches (mm): 15" x 30 (384 x 760)

Chamber Volume: 22Gal - 85 Liters

Overall Dimensions DxWxH - inches (mm): 34.5" x 26" x 20.6" - (875 x665 x 525)

Standard Cassettes Capacity:15 full

Tray Dimensions - inches (mm): 26.6" x 11.3" x 1" - (675 x 286 x 25), 26.6" x 13.8" x 1" - (675 x 286 x 25)

No. of Trays: 2

Standard Unwrapped Cycle Time - Cold: 32 Min. Hot: 23 Min.,/li>

Voltage (V) - Frequency (Hz): 230 - 50/60

Power (W) - Current (A): 3000W - 13.6A

Autoclave Weight - Lbs (kg): 213 (97)