



The Tuttnauer EZ9 Plus is a completely automatic autoclave. It was designed to have a modern layout with an easy to read Vibrant blue LED Viewing panel. The EZ9 Plus includes the newly designed EZglide self-locking door, which is ergonomically friendly and can be opened and closed effortlessly. The Autoclave is one of the largest mid-sized models on the market with a 9 inch by 19.8 inch chamber. There are now 3 ways to check on the cycle information, these options include a printer, a USB thumb drive and an ethernet cord that connects to your local network. Unlike other Autoclaves the EX-Plus Series has more tray surface to fit more instruments in the sterilizer.

Features

- The newly designed, ergonomically friendly EZGlide door features a sleek modern look and effortless opening and closing
- Innovative EZView multi-color display offers real time cycle information that is easy to read, even from a distance; the EZView displays cycle parameters, temperature, and pressure in the chamber; during a cycle each stage is displayed with an active color wheel and text
- Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility and ensure efficient drying; a high efficiency air pump circulates filtered air through the chamber for fast drying. 0.2 μ HEPA filter (replaceable)
- Comprehensive 2 year parts and labor warranty and 10 year chamber warranty
- 11 Inch Diameter X 19.8 Inch deep chamber and flexible loading options; the exceptionally deep 19.8 Inch chamber works in a variety of loading options to suit any office
- Use the 5 full sized trays for packs/pouches
- Use the tray rack (included) to accommodate 4 full plus 4 half sized cassettes no additional cassette rack to buy!
- Use the optional Tuttnauer pouch rack (model AR910) to encourage proper loading; one pouch rack is included with each EZ11Plus



Specifications

Dimensions

Height: 15.1 in (384 mm)
Width: 19.1 in (485mm)
Depth: 24.8 in (630 mm)
Weight: 104 lbs (47 kg)

Chamber Dimensions

9 in x 19.8 in (230 mm x 504 mm)
Chamber Volume: 5.2 gal - 19.8 liters
Tray Dimensions
Height: 0.6 in (16 mm)
Width: 6.7 in (170 mm),/li>
Depth: 16.3 in (415 mm)
No. of Trays: 3 trays

Power

Current: 1400 w - 12A
Voltage - Frequency: 120V - 50/60 HZ

