

Tuttnauer 2340M Autoclave



Steam Autoclave by Tuttnauer, the Tuttnauer 2340M Autoclave features a 9" Chamber. This Autoclave is designed to complement any healthcare facility. The Tuttnauer 2340's bright new easy to use panel was designed with the operator in mind and it will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

Features:

- Automatic shut off at the end of both the sterilization and dry cycles.
- A long life electro polished chamber and door.
- Double safety locking device prevents door from opening while chamber is pressurized.
- Overheat cutoff prevents burnout of heating elements in case of insufficient water in the chamber.
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir.
- Dual safety thermostat to protect against overheating.
- International certifications.



Typical manufacturer's specifications and images

Tuttnauer 2340M Autoclave

Specifications:

Chamber Dimensions:	9" x 18.5"
Volume:	5 gal
Overall Dimensions:	D x W x H: 21.5" x 20" x 14.4"
Standard Cassettes:	2 full/2 half
Tray Dimensions	D x W x H: 6.7" x 16.3" x 0.8"
No. of Trays:	3

Standard Unwrapped Cycle Time

- Cold: 30 mins
- Hot: 20 mins

Voltage (V):	120V
Freq. (Hz.):	50/60
Power (W):	1400W
Current (A):	12A



Soma Technology Inc.
166 Highland Park Drive
Bloomfield, CT 06002 • USA
1.800.438.7662
www.SomaTechnology.com

Note: The technical data given in this publication is for general information and are subject to change without notice. Actual configuration on the unit may vary. Contact our sales representatives for a complete list of details.