



The **2540EKA** autoclave significantly increase your productivity, providing extra fast and efficient drying cycles. This unit has the added benefit of a high efficiency air pump which allows closed door active drying. The 2540EKA is built for improved sterilization with the ability to dry packs and pouches.

The Tuttnauer 2540EKA Autoclave has a chamber that measures 10" x 19", with support for up to 4 trays each measuring 16.3" x 6.7". With the 4 trays available, users can quickly and efficiently sterilize equipment, with settings for sterilizing liquids, soil, waste products, instruments, and glassware. This particular model is equipped with a high efficiency air pump that offers an enhanced drying cycle.

The front panel is designed to be user-friendly, with programmable controls and a safety door that automatically seals when pressurized. The entire unit is made from a strong stainless steel construction. With this flexibility, the sterilizer is ideal for busy clinics, surgery centers, laboratories, dental offices, tattoo shops and other areas with a high-volume of patients or equipment use.

Features

- Enhanced drying cycle with a high efficiency air pump.
- Includes a replaceable 0.2 μ m Hepa air filter.
- Front panel programmable controls.
- Safety door automatically seals when pressurized and cannot open until depressurized. Closed door drying.
- Hepa filter 0.2 air filter (replaceable)
- Circulates hot air after the sterilization cycle for fast drying.



Specifications

Dimensions

Height: 14.4 in (36.6 cm)
Width: 20 in (51 cm)
Depth: 21.5 in (41 cm)
Weight: 81.6 lbs (37 kg)

Chamber Dimensions

Diameter: 10 in (25 cm)
Depth: 19in (25 cm)

Trays

Number of trays: 4
Height:14.4 in (36.6 cm)
Lenght: 21.5 in (55cm)
Width: 14.4 in (51 cm)

Cassette

Standard Cassettes Capacity: 3 full / 3 half.

Electrical

Voltage: 230 V.
Frequency: 50 / 60 Hz.
Current: 9.6 A.
Power: 2200 W.

Door Swing

Right Side Hinge