



The **Winco 6570** Drop Arm CareCliner has arms that can be lowered flush with the seat to aid in patient transfer. It has non-porous blow molded fold down side trays on both sides. These bariatric recliners have a 450 lb weight capacity. Each chair in the CareCliner Series features a standard Trendelenburg release with handle accessible from either side. It also has a built-in ergonomic headrest, adjustable elastic bands on the back frame that provides exceptional comfort and an elastic open seat base. All of these features combine to create the ultimate in comfort, durability, and value. For use in Dialysis, Oncology, Acute Care & other clinical environments.

Features

- Standard headrest cover
- Padded drop-down armrests
- dual fold down side-tables
- Easy-to-clean side panels
- Trendelenburg Positioning



Specifications

General

Weight: 111 lbs. (50.35 kg.)

Height: 49 in. (124.46 cm.)

Width with Tray Down: 37.75 in. (95.89 cm.)

Upright Depth: 43 in. (109.22 cm.)

Seat Height: 21 in. (53.34 cm.)

Seat Width: 25 in. (63.50 cm.)

Seat Depth: 20 in. (50.80 cm.)

Seat Back Height: 31.5 in. (80.01 cm.)

Armrest Height: 27 in. (68.58 cm.)

Weight Capacity: 450 lbs. (204.12 kg.)

Footprint Recline Length: 76.5 in. (194.31 cm.)

