



OPMI® 1 FC symbolizes quality, precision and reliability. By helping define the field of microsurgery, the system has built its own legacy.

The manual system is easy to use and delivers high fidelity images with the legendary Carl Zeiss optics. The precise fine focus functionality and the brightness adjustment knob allow the surgeon to quickly and easily attain the ideal focus and illumination levels - in routine examinations and surgical procedures.

After completing the procedure, the system can be conveniently moved away from the operating field and automatically shut off by raising the surgical microscope arm upward into the stand-by position.

#### Features

- Mobile OP-Microscope
- Oculare 10x22B
- adjustable view
- Zoom: 0,4-2,5
- Max. working distance 2
- Magnification 4,3x until 2
- Halogen light source 100
- Stand S2



Typical Manufacturer's Image