

Features

-
- A red medical cart with a white top tray and base. It features six drawers on the left side, each with a silver handle. On the right side, there is a vertical panel with a red heart-shaped button and a yellow latch. The cart is mounted on four wheels, with the front wheels being swivel casters.

Specifications

Dimensions

Height: 44 in (111.8 cm)

Width: 26.4 in (67.1 cm)

Depth: 19.2 in (48.7 cm)

Drawers

5

Drawer Heights

3 Drawers: 4.3 in (10.9 cm)

2 Drawers: 5.3 in (13.5 cm)

1 Drawer: 7.8 in (19.8 cm)

Features

Extension table

CPR board

IV Pole

Oxygen tank holder

External power strip

ABS Drawer dividers

Central or seals lock

Noise resistant castros with brakes

Comfortable push handle

Easy glide drawers

