

ROSS FLEXIFLO PATROL



ABOUT

The Ross Flexiflo Patrol incorporates the use of a specially-designed administration set which automatically clamps off the flow of fluid when removed from the pump, preventing free-flow. The Flexaflo Patrol features include a dose limit, accumulated volume counter, programmable hold time, lock-out, and memory which retains the rate, dose, and accumulated volume. The Patrol is recommended for pediatric patients. Rotary peristaltic operation.



FEATURES

Unique Safe-T-Valve feature on feeding set designed to prevent accidental product free flow if the set is improperly loaded or becomes dislodged from the pump

Free-flow alarm alerts caregiver to free-flow condition

Lock-out feature designed to prevent unauthorized changes

Special holder for safe, convenient storage of feeding set cover when not in use

Flow rates of 1 to 300 mL/hr in 1-mL increments

Dose feature 1 to 9999 mL in 1-mL increments

Memory retention for RATE, DOSE, and VOLUME FED

Unique programmable HOLD of 1 to 90 minutes

Easy-to-use dial control and display panel

SPECIFICATIONS

DIMENSIONS

Height: 8.5 in

Width: 6.5 in

Depth: 4.8 in

Weight: 6.6 lbs

POWER

Voltage: 100 to 130 VAC, 60 Hz, 1-Phase

Power: 15 watts, average

Fuse: 200 mA, time-delay type

Line Cord: Hospital-grade, non-detachable (10 feet)

Leakage: Less than 100 microamps

SPECIFICATIONS



GENERAL

Rate Range: 1-300 ml/hr in 1 ml increments

Dose Limit: 1-9999 mls. In 1 ml increments

Accuracy: $\pm 10\%$ or ± 0.5 ml, whichever is greater

Occlusion: 15-23 PSI

Battery Life: Approx. 3 hrs. at 125 ml/hr

Recharge: Max. or 4.5 hrs. while the pump is off or 12 hours during use

Alarms: Free flow, No flow, Dose fed, Lo Batt, Select Rate, Select run, Internal Malfunction (E1-E9)

Memory: Auto retention of rate, dose, and accumulated volume.