

MEDIVATORS SCOPE BUDDY



ABOUT

The **Medivators Scope Buddy** is an endoscope flushing aid that automates the channel-cleaning step of endoscope reprocessing. Instead of manually syringing each channel, staff press one button and the unit delivers a steady, timed flow of fluid through every channel until the cycle finishes. It works with all flexible endoscopes, including ultrasound scopes and ones with elevator wire channels, and it can flush more than one scope at a time. A full cycle takes about 35 seconds, and the unit signals with a sound and a light when it's done. Because it takes over the syringing that's normally done by hand, it also cuts down on the repetitive strain that comes with doing it manually.

FEATURES

Hands-Free Flushing: Replaces manual syringing of channels.

One-Touch Operation: A single push-button starts the cycle, no manual timing needed.

35-Second Cycle: Thoroughly flushes most endoscopes in about 35 seconds.

Broad Compatibility: Works with all flexible endoscopes, including ultrasound scopes and those with elevator wire channels.

DIMENSIONS

Height: 9.5 in (24.1 cm)

Width: 6 in (15.2 cm)

Depth: 8 in (20.3 cm)

Weight: 10 lb (4.5 kg)



OPERATION

Fluid Delivery Rate: 1 liter/minute

Flush Time: 0 - 99 minutes

Timer Control: Digital display with touchpad controls

Cycle Completion Indicator: Audible and visual alerts

System Pressure: Built-in safety bypass