

ECOLAB ORS-1075HS



ABOUT

The **Ecolab ORS-1075HS** hush slush machine is a workhorse in the OR, churning out ultra-smooth surgical slush to keep tissues and organs at the right temperature during procedures. It automatically mixes the slush so it stays perfectly consistent, no sharp ice chunks that could harm delicate tissue. The Disc-Draper tech makes setup and cleanup simple, so your team spends less time prepping and more time focused on surgery. With a 7.5-liter capacity and optional fluid warming, it's flexible enough to handle whatever the case demands. Dependable, efficient, and built for real-world OR challenges.

FEATURES

- **Lush Technology:** Produces soft, velvety slush ideal for tissue and organ cooling during surgery.
- **Automatic Stirring:** Prevents the formation of sharp ice edges, ensuring patient safety and consistent slush quality.
- **7.5-Liter Basin Capacity:** Provides ample volume to meet the needs of demanding surgical procedures.
- **Disc-Draper Technology:** Streamlines preparation and cleanup for enhanced operating room efficiency.
- **Optional Fluid Warming System:** Adds versatility by maintaining or warming fluids as needed.



SPECIFICATIONS



DIMENSIONS & WEIGHT

Main Unit

- Height: 36.0 in (91.4 cm)
- Width: 18.0 in (45.7 cm)
- Depth: 18.0 in (45.7 cm)
- Weight: 98 lbs (44.5 kg)

Rear Storage Pod

- Height: 9.8 in (24.8 cm)
- Width: 17.8 in (45.2 cm)
- Depth: 7.0 in (17.8 cm)

FLUID CAPACITY

Basic Reservoir Capacity: 7.5 Liters (Round)

Operating Fluid Level Range: 1 to 4 Liters