

STRYKER MX-PRO



ABOUT

The Stryker MX-PRO is lightweight, durable, and easy to maneuver. This manual cot has a weight capacity of 650 lbs while only weighing 81 lbs. It combines strength and durability with ease of handling. The one-hand operation for the backrest, side rail, breakaway head section, and safety bar allows the user to focus on life-saving tasks. The PRO Stretcher features a foot-end lift and grip design, along with duplicate foot-end controls, which provide ergonomic support for operators of varying sizes.



FEATURES

- Scientifically optimized lift bar design
- Positive action height adjustment
- One-hand release breakaway head section with safety bar
- Oversized wheels maximize mobility on all types of terrain
- The wheel lock option helps keep the unoccupied cot from rolling
- One-four-point shoulder harness restraint
- Two lap belt restraints
- One-hand release, fold down side rails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- Fold-down side rails enable caregiver-patient access even in the smallest spaces.
- X-frame that features 7 height positions



SPECIFICATIONS



HEIGHT

- Position 1: 13.5" (34cm)
- Position 2: 21.0" (53cm)
- Position 3: 25.5" (65cm)
- Position 4: 29.0" (74cm)
- Position 5: 32.0" (81cm)
- Position 6: 35.0" (89cm)
- Position 7: 37.5" (95cm)

LENGTH

- Standard: 80.5" (205 cm)
- Minimum: 62.0" (158 cm)

WIDTH

- 23" (58 cm)

WEIGHT

- 181 lbs (38kg)

WHEELS

- Diameter: 6" (15 cm)
- Width: 2" (5cm)

ARTICULATION

- Backrest: 2 - 73°
- Shock Position: +14°

MAXIMUM WEIGHT CAPACITY

- 650 lbs (295 kg)

MINIMUM OPERATORS REQUIRED

- Occupied Cot: 2
- Unoccupied Cot: 1

RECOMMENDED LOADING HEIGHT

- Up to 32" (81cm)