

ECOLAB ORS - 1058HS



ABOUT

The Ecolab ORS-1058HS hush slush machine delivers a velvet-soft slush during procedures that require it. The ORS-1058HS utilizes proprietary Auto-Stir technology to stir the solution while cooling automatically. This eliminates concerns with conventional slush that may have sharp edges, which can cause trauma to delicate tissues. Check out the ORS-1075HS for a larger basin capacity.

FEATURES

Auto-Stir technology: Automatically stirs and cools the solution to create consistently soft surgical slush, eliminating the sharp ice edges common in conventional systems.

Disc-Drape system: Snaps easily onto the Auto-Stir device, maintaining optimal sterile coverage and ensuring a safe surgical environment.

Locking caster wheels: Provide smooth mobility for easy transportation while locking mechanisms minimize accidental movement during procedures.

Quiet operation: Reduces noise in the operating room, allowing staff to maintain focus and clear communication.

Optional ORS Fluid Warming System: Offers flexibility for facilities needing both slush creation and fluid warming capabilities in one integrated solution.



SPECIFICATIONS

DIMENSIONS

Height: 36" (91.44 cm)

Width: 18" (45.72 cm)

Depth: 18" (45.72 cm)

BASIN CAPACITY

5.8 Liters (Rectangle)

FLUID LEVEL

1 to 3 Liters