

AXIA PRO-SERIES PROCESSING SINK



ABOUT

The **Axia Processing Sink** offers a variety of models to fit your facilities needs. It features a fully welded construction with a reinforced top so you know you're getting a quality stainless steel surgical sink. It can be configured to multiple different layouts depending on individual requirements. This Axia Processing Sink can be purchased in single, double, and triple washing stations. The standard features of this scrub sink are listed below.

FEATURES

- 12" Angled Backsplash
- 37" Standard Height for Fixed Height Sinks
- Adjustable Feet
- Basin Size: 24"W X 20"L X 10"D
- Etched Gallon Markings
- Flat Working Surface
- Fully Welded Construction with Reinforced Tops
- Marine Edge
- Open Pipe Base
- Seismic Brackets
- Single, Dual, or Triple Basins



SPECIFICATIONS



OPTIONAL MODEL FEATURES

Ergonomically friendly electronic adjustable height hydraulic system with four position programming (electronic adjustable height range: 35"- 42")

Bolt-on Side Splash

Copper Drain System (Not Compatible with Deionized Water)

Hand Lever Drains

Larger Basin: 30"W X 20"L X 10"D

OPTIONAL SINK ACCESSORIES

Backsplash Mounted Accessories

- Air Gun Kit (with 5 Tips)
- Deionized (DI) Water Faucet
- Eye Wash
- Gooseneck Mixing Valve Faucet
- Pre-Rinse Faucet
- Water Gun Kit (with 8 Tips)

Basin Covers

Leg Mounted Power Strip

Lower Shelf/Shelves

Wrist Rest

OPTIONAL PEGBOARD BACKFRAME

ACCESSORIES

Bin Rail

- Plastic Storage Bins – 2 Sizes Available

Brush Holder

Digital Basin Temperature Indicator

Dosing Soap Pump

Flushing System

Gallon Jug Holder

Glove Dispenser

Keyboard Arm

LED Magnifying Inspection Lamp

LED Soak Timer

Monitor Mount Arm

Overhead LED Light

Overhead Light Deflector

Pegboard Hooks

Pegboard Shelf

Spray Bottle Holder

Waterproof LED Task Light with Magnifying Lens

Work Flow Signs

STANDARD MODELS AND DIMENSIONS

AX307237-PS – 30" D x 72" W x 37" H, Single Basin

AX309637-PS – 30" D x 96" W x 37" H, Dual Basin

AX3012037-PS – 30" D x 12" W x 37" H, Triple Basin

ADJUSTABLE HEIGHT MODELS AND DIMENSIONS

AX307235-PSE – 30" D x 72" W x 35" H, Single Basin

AX309635-PSE – 30" D x 96" W x 35" H, Dual Basin

AX3012035-PSE – 30" D x 120" W x 35" H, Triple Basin