

ZIMMER INTELLICART SYSTEM



ABOUT

The Zimmer IntelliCart waste management system sets the industry standard with its impressive 34-liter fluid capacity, ultra-quiet vacuum pump, obstruction-free suction manifolds, and portable smoke evacuation capability. Experience the ease and convenience of effectively managing surgical fluid and smoke waste hazards with the IntelliCart..

FEATURES

- Redesigned Vacuum Pump that is extra quiet with HEPA air filtration, increased flow capabilities, and odor control.
- Dual 17-liter lighted reservoirs with ConSeal™ Tint Technology selectively hide surgical fluids.
- Clog-free manifold with integrated high-flow filter protects against vacuum system moisture ingress.
- Simple and interactive, high-definition touchscreen display.
- Reusable reservoirs and recyclable manifolds reduce OR red bag waste by up to 70%

DIMENSIONS

Height: 55 in (139.7 cm)

Width: 23.5 in (59.6 cm)

Depth: 23 inches (58.4 cm)

DISPLAY

LCD, Backlit touch screen display



SPECIFICATIONS



SUCTION/EXHAUST

Closed suction system with variable suction levels based on user selection, providing regulated vacuum from 0-600 mmHg using AC-powered twin piston vacuum pump passing through single-use manifolds and into canisters..

IV POLE

Carts available with powered or manual IV poles, capable of 12,000mL maximum fluid capacity

SAFETY ALERT

Audible and visual display for early warning of overflow protection system and vacuum start-up alert

CANISTER EXTERIOR GRADUATIONS

Graduations marked on the canister exterior provide an approximate visual indication of contents in 100mL increments..