

AXIA FH-2000



ABOUT

The **Axia FH-2000** functional headwall is a versatile and efficient solution for patient care settings, designed for easy installation and maintenance. Its customizable configuration allows healthcare facilities to adapt it to various environments, including patient rooms, ICU, NICU, PACU, trauma, and L-D-R. The Axia FH-2000 comes pre-piped and wired to streamline setup and features essential utilities, including a nurse call system, multiple electrical outlets, and DISS gas outlets for oxygen, medical air, and vacuum. Its compact design optimizes space while maintaining functionality. The surface-mounted structure ensures durability, providing reliable support for medical equipment and supplies.



FEATURES

Easy Installation and Maintenance: Pre-wired and pre-piped for straightforward setup and efficient upkeep.

Integrated Nurse Call System: Enhances communication and care responsiveness.

Multiple Utility Outlets: Includes 2 ivory and 1 red electrical outlets, along with DISS outlets for oxygen, medical air, instrument air, and vacuum. It also includes a slide bracket.

Space-Efficient Design: Surface-mounted with a length of 48 in (122 cm), height of 8.7 in (22 cm), and depth of 2.6 in (6.6 cm).

DIMENSIONS

Height: 8.7 in (22 cm)

Width: 48 in (122 cm)

Depth: 2.6 in (6.6 cm)