

# STRYKER POWER PRO IT



## ABOUT

The Stryker Power PRO-IT is a powered incubator transport stretcher designed to enhance caregiver safety and efficiency while ensuring patient stability during transport. With its innovative battery-powered hydraulic system, the Power PRO-IT minimizes manual lifting, reducing the risk of caregiver injuries and streamlining patient handling in emergency and critical care environments.

## FEATURES

- **Battery-Powered Hydraulic Lift System:** Effortlessly raises and lowers the stretcher with the touch of a button.
- **Automatic High-Speed Retract:** Saves time during loading and unloading.
- **Manual Backup System:** Provides confidence and reliability in all powered functions.
- **Adjustable Load Height:** Adapts to different ambulance configurations with load heights ranging from 26 to 36 inches (66–91 cm).
- **Durable Lightweight Design:** Constructed with a rugged aluminum frame and oversized wheels for stability and mobility.
- **Power-LOAD Compatibility:** Fully compatible with the Power-LOAD cot fastening system for enhanced safety and ease of operation.
- **Ergonomic Grip and Control System:** Scientifically optimized lift bar and control design for improved operator comfort.
- **Safety Features:** Includes a safety bar, tie-down brackets, and retractable head section for secure and versatile usage.



# SPECIFICATIONS



## DIMENSIONS

Height (Infinite Positioning):

Highest Position: 41.5 in (105 cm)

Lowest Position: 14 in (36 cm)

Length:

Standard: 63 in (160 cm)

Maximum: 81 in (206 cm)

Width: 23 in (58 cm)

Weight: 134 lbs (61 kg)

## WEIGHT CAPACITY

Maximum Capacity: 700 lbs (318 kg)

Unassisted Lift Capacity: 500 lbs (227 kg)

## WHEELS

Diameter: 6 in (15 cm)

Width: 2 in (5 cm)

## POWER SYSTEM

**Battery Type:** SMRT power system with 12V DC, 120V AC, and 240V AC options

**Compatibility:** Power-LOAD Model 6390, Floor Mount Models 6370 or 6377, Wall Mount Model 6371

## SAFETY AND CERTIFICATIONS

**Certification:** Conforms to IEC-60601-1 and KKK-A-1822 standards for safety.

### Warranty

- Two years on parts and labor (optional extensions available).
- Lifetime warranty on welds.