

# ETHICON MEGADYNE MINI VAC



## ABOUT

---

The **Ethicon MEGADYNE MINI VAC** is a portable smoke evacuation system designed for small procedure rooms and facilities that require effective plume management without the footprint of a larger unit. Despite its compact size, it delivers consistent performance, helping maintain a clear surgical field and safer environment for staff and patients. Its quiet operation and straightforward setup make it an efficient choice for practices seeking a simple yet dependable smoke evacuation solution.

## FEATURES

---

**Portable Size:** Lightweight and compact, making it easy to move between rooms or procedures.

**Reliable Smoke Capture:** Designed to effectively remove surgical smoke and maintain clear visibility during procedures.

**Low Noise Design:** Operates quietly to reduce disruption in office or clinical settings.

**Easy to Use:** Quick setup with intuitive connections for efficient operation.

**Affordable Option:** Provides effective smoke evacuation performance while minimizing costs for facilities.

## DIMENSIONS

---

Height: 7.9 in (20 cm)

Width: 7.9 in (20 cm)

Depth: 10.8 in (27.5 cm)

Weight: 6.8 lbs (3.1 kg)



## PERFORMANCE

---

Maximum Flow Rate: Up to 35 LPM

Designed for localized smoke evacuation during minor surgical procedures

Quiet operation (< 55 dBA)

## FILTRATION SYSTEM

---

ULPA filter with 99.999% efficiency at 0.12 microns

Effective removal of surgical plume, particulates, and odor

Filter life: approximately 6 hours of use

## POWER REQUIREMENTS

---

100–120 V, 50/60 Hz, 2.5 A

Standard hospital-grade power connection