

VERATHON GLIDESCOPE CORE 15 / 10



ABOUT

The **Verathon Glidescope Core 15** and **10** are versatile airway visualization systems designed for both video laryngoscopy and bronchoscopy. Available with a 15-inch or 10-inch Full HD touchscreen monitor, it delivers clear, detailed images to support confident airway management in both routine and challenging situations. The system's **Dual View** capability allows simultaneous visualization from a Glidescope video laryngoscope and a bronchoscope, helping clinicians act quickly and precisely. With its ergonomic workstation, organized storage, and advanced imaging features, the Glidescope Core streamlines workflow while enhancing safety and efficiency.

FEATURES

Beautiful visualization with Full HD display, Dynamic Light Control, and anti-fog technology for consistently clear images.

Dual View capability for simultaneous video laryngoscopy and bronchoscopy to assist in complex airway procedures.

Mobile workstation with articulating monitor arm, prep tray, storage bins, and organized tool access.

Compatible with a wide range of single-use and reusable Glidescope video laryngoscopes and BFlex bronchoscopes in multiple sizes.

Infection control friendly design with integrated disposal bin and protective scope sleeve.

MONITOR SIZES

Core 15: 15" Full HD touchscreen

Core 10: 10" Full HD touchscreen



VIEWING MODES

Dual View (Core 15) – side-by-side live image display

Picture-in-Picture (Core 10) – inset image display during dual procedure

INPUT/OUTPUT

Dual scope input ports for seamless switching or concurrent viewing

HDMI out for external video display support