

BAXTER HILLROM TOTALCARE P1900



ABOUT

The Baxter Hillrom **TotalCare P1900** Bed System is a comprehensive, modular acute care bed designed for maximum versatility across care environments, including ICU, critical care, step-down, medical/surgical, and emergency departments. It integrates advanced mechanical, hydraulic, pneumatic, and electronic systems to provide enhanced patient safety, caregiver efficiency, and therapeutic functionality.

FEATURES

Modular Design: Configurable with various optional modules to suit specific patient and facility needs.

Advanced Positioning: Offers a wide range of positions including Chair, Chair Egress, Recliner, Trendelenburg, Reverse Trendelenburg, and CPR.

Patient Safety Systems: Includes bed exit detection, nurse call integration, siderail controls with lockout features, and a comprehensive alarm system.

Caregiver Efficiency: Features Point-of-Care siderail controls, HandsFree emergency releases, and optional powered transport (IntelliDrive).

Pressure Management: Compatible with advanced air surfaces for pressure injury prevention.

SAFE WORKING LOAD

500 lb (227 kg) (patient and equipment)

HEIGHT RANGE

Low Position: Between stepped sides, 14.5 in – 15.5 in (36.8 cm – 39.4 cm)

High Position: Between stepped sides, 33.5 in – 34.5 in (85.9 cm – 87.6 cm)



OVERALL WIDTH

Siderails Raised: 40.5 in (102.9 cm)

Siderails Lowered: 37 in (94 cm)

Distance Between Casters: 25.5 in (64.8 cm)

OVERALL LENGTH

Overall Length: 91.5 in (232.4 cm)

With Roller Bumpers Installed: 92.5 in (235 cm)

With Transport Shelf: 93.5 in (237.5 cm)

FOOT SECTION DIMENSIONS

Width: 28.88 – 29.13 in (73.36 – 73.99 cm)

Length: 26.75 – 27.25 in (67.36 – 73.99 cm)

SPECIFICATIONS



MINIMUM UNDER BED CLEARANCE

- 4.75 in (12.07 cm)
- 4.25 in (10.80 cm) ("B" model beds and newer)
- 1.25 in (3.18 cm) (IntelliDrive Transport System Beds)

ARTICULATING DECK LENGTH

Length: 81.75 – 82.25 in (207.65 – 208.92 cm)

MOVEMENTS

- Maximum Head Elevation:** 75°
- Maximum Knee Elevation:** 20°
- Maximum Foot Elevation:** 85° down
- Maximum Foot Retract:** 12 in (30 cm)
- Maximum Trendelenburg:** 15°
- Maximum Reverse Trendelenburg:** 15°
- Maximum Emergency Trendelenburg:** 20°

BED WEIGHT

- With Minimum Options:** 525 lb (238 kg)
- With Maximum Options:** 585 lb (265 kg)

CASTERS

Diameter: 5 in (13 cm)

SLEEP SURFACE LENGTH

- 3 in (211 cm) Total
- 56 in (142 cm) Primary Section
- 27 in (69 cm) Foot Section

SLEEP SURFACE WIDTH

- 35 in (89 cm) Primary Section
- 31 in (69 cm) Foot Section

SCALE

Maximum Scale Weight: 400 lb (181 kg)