

SMITH AND NEPHEW T-MAX BEACH CHAIR



ABOUT

The Smith and Nephew **T-MAX Beach Chair** positioner is designed to provide optimal support and accessibility for both arthroscopic and open shoulder procedures. It enables the entire surgical team to position patients securely and efficiently in the beach chair configuration. The system features a dual ball and socket head rest for superior anesthesia access and patient head stability. Additionally, its sliding shoulder pads and lateral supports enhance exposure and support during operations. The T-MAX is compatible with a variety of operating tables, ensuring versatility across different surgical environments.

FEATURES

Hydraulic assist mechanism for smooth raising and lowering of the patient

Sliding shoulder pad and lateral adjustment to optimize access to the surgical site

Integrated intubation tree and a wide range of head rest movement for anesthesia management

Adjustable lateral supports and knee bolster to ensure patient stability and prevent sliding

A dedicated spider rail for compatibility with the SPIDER2 limb positioner

SAFE WORKING LOAD

Maximum weight limit for O.R. table must not be exceeded. Most older tables have a 300 lb (136 kg) limit. The T-MAX will accommodate patients up to 350 lbs (159 kg).

