

# BAXTER HILLROM TOTALCARE BARIATRIC +



## ABOUT

The **Baxter Hillrom TotalCare Bariatric Plus Therapy System** is a heavy-duty critical care bed built for the needs of larger patients. It supports up to 500 lbs (227 kg) and combines mobility, pressure relief, and positioning tools that help reduce caregiver strain and maintain patient comfort. The bed's surface options help protect skin health, while powered mobility features make transport easier for hospital staff.

## FEATURES

- **MicroClimate Management:** Uses Low Air Loss surfaces to help maintain dry skin conditions and offers options such as turn assist, percussion and vibration, and continuous lateral rotation therapy.
- **Motorized Transport:** The IntelliDrive XL system moves patients through the facility with minimal effort from caregivers.
- **Chair Position & Seat Deflate:** Converts smoothly from a flat surface to a seated position and includes a seat deflate mode to assist with transfers or repositioning.
- **Built-In Weighing Scale:** Provides weight readings within 1% accuracy (up to 400 lbs) and handles patients up to 500 lbs.
- **Adjustable Length:** Extends from 72 to 84 inches to fit different patient heights.
- **Safety Features:** Offers bed exit alerts, head angle alarms, and clear visual indicators to support safe patient handling.
- **Siderail Controls:** Simple controls on the siderails allow quick adjustments without leaving the bedside.



# SPECIFICATIONS



## DIMENSIONS

### Overall Bed

- **Maximum Overall Width:** 44.5 in (113 cm)
- **Maximum Overall Length:** 81.5 in – 93.5 in (206.8 – 237.3 cm)
- **Height (Sleep Deck to Floor):** Low: 19 in (48.2 cm), High: 38 in (94 cm)
- **Weight Capacity:** 500 lbs (227 kg)

### Pressure Redistribution Foam

- **Length:** 72-84 in (182.7-213.4 cm)
- **Width:** 40 in (102 cm)
- **Height:** 7 in (17.8 cm)

### Microclimate Management:

- **Length:** 72-84 in (182.7-213.4 cm)
- **Width:** 40 in (102 cm)
- **Height:** 11 in (27.9 cm)

## POWER

**Power Requirements:** 120v, 60 cycle AC

**Fuse, Maximum Current:** 15A, 9.9A

**Safety Listing:** UL60601-1

## BED ANGLES

**Head Section:** 0° – 75°

**Knee Section:** 0° – 20°

**Max Trendelenburg:** 15°

**Reverse Trendelenburg:** 20°

## CASTER & BRAKING SYSTEM

**Caster Diameter:** 5 in (12.7 cm)

**Four Wheel Braking:** Yes

**Steer and Lock Both Rolling and Swivelling:** Yes

## SCALE SYSTEM

**Accuracy:** 1% of patient weight (up to 400 lbs)

**Scale Weight Capacity:** 500 lbs (227 kg)