

PEDIGO P-2120



ABOUT

The **Pedigo P-2120** is a dual-basket fluid warming cabinet built to heat irrigation and injection solutions to safe clinical temperatures. It uses an electronic control panel with dedicated settings for different fluid types, helping maintain proper temperature control during use. The cabinet has a stainless steel body with a clear door for quick inventory checks. A built-in fan circulates air evenly, while the epoxy-coated baskets and sturdy construction make it a dependable option for everyday clinical work.

FEATURES

- **Dual-Basket Capacity:** Two epoxy-coated wire baskets hold up to 9 one-liter bottles or 15 one-liter bags each, for a total of 18 bottles or 30 bags.
- **Dual Warming Modes:** Separate temperature ranges for irrigation fluids (37°–66°C / 98°–150°F) and injection fluids (37°–40°C / 98°–104°F).
- **Forced-Air Circulation:** An internal fan keeps air moving throughout the chamber to maintain even heat.
- **Digital Controls:** Simple LED control panel with mode buttons and visual or audible alerts for over-temperature protection.
- **Easy Transport:** Four 5-inch casters, two with brakes, and a perimeter bumper for safe movement around the facility.
- **Built for Safety:** Classified as Safety Class I, Hospital Grade equipment, with an independent thermal shut-off system and compliance with medical standards.



SPECIFICATIONS



DIMENSIONS

Width: 28.53 in (724.5 mm)

Height: 45 in (1148 mm)

Depth: 24.13 in (612.7 mm)

POWER

Voltage: 120 VAC

Frequency: 50/60 Hz, 1 Phase

Power Consumption: 0.78 kW, 6.4 Amps

Plug Type: NEMA 5-15P

Mode of Operation: Continuous

ENVIRONMENTAL CONDITIONS

Operating Temperature: 15°C to 32°C (60°F to 90°F)

Storage Temperature: -40° to +70°C (-40° to +159°F)

Storage Relative Humidity: 10% to 95%, non-condensing

CAPACITY & PERFORMANCE

Basket Capacity: 2 Baskets, each holding 9 x 1L bottles or 15 x 1L bags

Temperature Accuracy: Within +0/-1.12°C (+0/-2°F) for 37°-40°C range; Within +0/-1.67°C (+0/-3°F) for 43°-66°C range