

PEDIGO P-2220



ABOUT

The **Pedigo P-2220** is a high-capacity **blanket warmer** designed to provide consistent, safe warmth for patient comfort in demanding clinical environments. Built with reliability in mind, it includes safety features like an independent warming shut-off system and power-fail recovery to ensure consistent performance in busy hospitals and surgical centers.

FEATURES

- **WarmSafe Heating:** Employs patented technology that heats specific zones as needed, monitored by multiple sensors for optimal efficiency and temperature uniformity.
- **Substantial Capacity:** Offers ample space for storing and warming a large volume of cotton blankets or towels in a single unit.
- **Advanced Control:** Features a four-digit LED display, adjustable temperature from 90°F to 160°F, and a control lockout function to prevent unauthorized changes.
- **Safety Systems:** Includes a separate warming shut-off system to prevent overheating and retains programming to resume normal operation after a power interruption.
- **Heated Center Shelf:** Includes a dedicated heated shelf to further optimize heat distribution and performance throughout the entire cabinet cavity.
- **Mobility Features:** Full-perimeter bumper and four 5-inch heavy-duty casters, two with locks.
- **Air Circulation:** Fans in each chamber move air to avoid hot or cold spots within the cabinet.



SPECIFICATIONS



DIMENSIONS

Height: 28.0 in (711 mm)
Width: 18.5 in (470 mm)
Depth: 23.0 in (584 mm)
Weight (Net): 91 lbs (41 kg)

POWER

Input Voltage: 120 VAC
Input Frequency: 60 Hz, 1 Phase
Power Consumption: 0.8 kW, 6.7 Amps
Safety Classification: Class I Equipment, No Applied Parts

ENVIRONMENTAL CONDITIONS

Operating Temperature: 60°F to 90°F (15°C to 32°C)
Operating Humidity: Above 20%, non-condensing
Storage Temperature: -40°F to 159°F (-40°C to +70°C)

CAPACITY

Internal Volume: 3.5 ft³