

# AXIA SB5



## ABOUT

The **Axia SB5** pain management table is a five-movement imaging table designed to support a wide range of pain management and minimally invasive procedures. Built with a carbon fiber tabletop and a compact H-style base, the SB5 provides excellent C-arm access while maintaining stability across all patient positions. Its generous surface area and high weight capacity allow clinicians to work comfortably with a broad patient population, including bariatric cases. Smooth powered movements and a low minimum height help simplify patient positioning and transfers..



## FEATURES

- **Five Powered Movements:** Includes height adjustment, lateral tilt, Trendelenburg, longitudinal travel, and lateral travel.
- **Carbon Fiber Imaging Top:** Radiolucent, diving-board style tabletop allows unobstructed C-arm access.
- **High Patient Capacity:** Rated to support patients up to 500 lb., regardless of table position.
- **Large Patient Surface:** An 84-inch by 24-inch tabletop accommodates taller patients and eliminates the need for arm boards.
- **Low Entry Height:** Designed to make patient boarding easier and safer for both staff and patients.
- **Smooth Mobility and Locking:** Full-swivel casters with reliable wheel locks allow easy repositioning within the procedure room.
- **Optional Cervical Top Configuration:** Available with an optional cervical top to support cervical spine and upper pain management procedures.



\*Optional Cervical Top

# SPECIFICATIONS



## TABLETOP DIMENSIONS

**Length:** 84 in (214.4 cm)  
**Width:** 24 in (61 cm)  
**Thickness:** 1.62 in (4.1 cm)

## RADIOLUCENT DIMENSIONS

**Length:** 56 in (142.2 cm)  
**Width:** 24 in (61cm)

## TABLE DIMENSIONS

**Length:** 86 in (218.4 cm)  
**Width:** 26 in (66 cm)

## WEIGHT

**Table Weight:** 840 lb. (381 kg)  
**Patient Capacity:** 500 lb. (227 kg)

## TABLE MOVEMENTS

**Height:** 33.5 to 43.5 in (85 to 110.5 cm)  
**Lateral Tilt:** ( + / - ) 17°  
**Trendelenburg:** 13°  
**Longitudinal Travel:** 10 in (+/- 5 ")  
**Lateral Travel:** 10 in (+/- 5 ")

## GENERAL

**Mattress Thickness:** 2 in (5 cm)  
**Tabletop / Mattress Attenuation:** 1.6mm Al. (+/- .2mm tolerance)  
**Tabletop Material:** Carbon fiber skin/foam core with integral facial cut-out  
**Casters:** 5 in casters with full swivel and wheel lock  
**Steel Base:** H-Shaped with epoxy finish  
**Hand Control:** Activates all motions  
**Patient Straps:** 2

## POWER

120V/60HZ or 230V/50HZ  
Protect product from moisture  
**Allowable Storage & Operating Humidity:** 0-92%  
**Allowable Storage Temperature:** 32°F to 104°F (0°C to 40°C)  
**Allowable Operating Temperature:** up to 104°F (40°C)