

# HILLROM YELLOFINS APEX



## ABOUT

The **Hillrom Yellofins Apex** stirrups are designed to provide stable, controlled lithotomy positioning while supporting both patient safety and surgical team efficiency. Built for demanding surgical environments, these reusable stirrups offer smooth positioning adjustments directly from the sterile field, helping reduce repositioning time and physical strain on staff. The design supports a wide range of procedures, including urology, gynecology, colorectal, and robotic-assisted surgeries. With an emphasis on infection control and ease of use, the Yellofins Apex stirrups integrate seamlessly into modern operating rooms. The result is consistent positioning support that helps teams focus on the procedure rather than equipment limitations.



## FEATURES

- **High Patient Weight Capacity:** Supports patients up to 500 lb
- **Sterile-Field Adjustability:** Allows lithotomy, abduction, and lateral adjustments to be made from the sterile field.
- **Dual-Bar Safety Design:** Engineered to help prevent medial leg drop.
- **Single-Point Boot Release:** Enables quick, controlled leg positioning adjustments with minimal effort.

## SPECIFICATIONS

- Patient weight capacity:** 500 lb
- Lithotomy adjustment (from sterile field):** Yes
- Lithotomy range (from sterile field):** +80° to -32°
- Abduction range (from sterile field):** +23° to -9°
- Maximum medial/lateral shift:** ±20°
- Product weight (single stirrup):** 18.5 lb (8.4 kg)