

3M RANGER PRESSURE INFUSOR



ABOUT

The **3M Bair Hugger Pressure Infusor Model 145** is designed to support rapid fluid delivery when time matters most. Commonly used in operating rooms, trauma settings, and other critical care environments, this unit applies consistent pressure to IV solution bags to help clinicians deliver fluids quickly and efficiently. The system works with standard IV fluid bags ranging from 250 mL to 1000 mL and is compatible with a wide range of solutions, including blood, saline, sterile water, and irrigation fluids. Built for adult and pediatric use, it provides dependable performance in high-demand clinical situations.

FEATURES

- **Rapid Infusion Support:** Designed to apply pressure to IV fluid bags when fast delivery of liquids is required.
- **Up to 300 mmHg Pressure:** Works with compatible administration sets capable of delivering fluids under pressures up to 300 mmHg.
- **Flexible Bag Compatibility:** Accepts IV solution bags from 250 mL to 1000 mL.
- **Broad Clinical Use:** Suitable for operating rooms, emergency trauma settings, and other acute care environments.
- **Compatible with Standard Fluids:** Supports use with blood products, saline, sterile water, and irrigation solutions that meet applicable medical standards.



DIMENSIONS

Height: 15.75 in (40 cm)

Width: 20 in (51 cm)

Depth: 7.75 in (20 cm)

Weight: 17 lb (7.7 kg)

PERFORMANCE

Operating Pressure: 300 ± 10 mmHg setpoint