

# PEDIGO P-2130



## ABOUT

The **Pedigo P-2130** fluid warmer is designed to safely warm and store irrigation or injection fluids in a consistent, controlled environment. Built with a stainless steel exterior and a single warming chamber, the unit is made to handle daily clinical use while maintaining even heat distribution throughout the cabinet. An internal circulation fan helps eliminate temperature stratification, ensuring fluids are warmed evenly and remain within tight temperature tolerances. An intuitive digital control system allows staff to easily set and monitor temperatures, with separate ranges available depending on the type of fluid being stored standards.

## FEATURES

**Three Basket Storage System:** Includes three epoxy-coated wire baskets, each designed to hold fluid bags or bottles, with capacity for up to 54 one-liter bottles or 72 one-liter bags total.

**Different Temperature Modes:** Programmable settings for irrigation (98°F to 150°F) and injection fluids (98°F to 104°F).

**Consistent Temperature Control:** Internal fan circulates air to maintain uniform chamber temperatures, with tight accuracy ranges across both operating modes.

**User-Friendly Controls:** Digital interface with LED display, simple adjustment buttons, and lockout feature to help prevent unintended changes.

**Built-In Safety Features:** Audible over-temperature alarm and independent thermal shut-off system help protect stored fluids and maintain safe operation.



## DIMENSIONS

**Height:** 40.8 in (103.6 cm)

**Width:** 22.3 in (56.6 cm)

**Depth:** 35 in (89 cm)

## BASKETS (X3)

**Height:** 7.7 in (19.5 cm)

**Width:** 17.1 in (43.4 cm)

**Depth:** 22.1 in (56.1 cm)