

ZIMMER DORNOCH TRANSPOSAL ULTRAFLEX FLUID CART

ABOUT

The **Zimmer Dornoch Transposal UltrafleX Fluid Cart** is a closed fluid waste management system designed to help streamline surgical fluid collection and disposal in the operating room. Built for efficiency, the UltrafleX Duo configuration utilizes dual reservoirs to handle high fluid volumes while minimizing staff exposure to potentially hazardous waste. Its mobile design allows for easy transport between procedure rooms, while the sealed system supports a cleaner and safer environment. The system is intended to simplify workflow by reducing manual handling and improving turnover times.

FEATURES

- **Dual Reservoir Capacity:** UltrafleX Duo configuration supports up to 33 liters of fluid collection.
- **Closed Waste System:** Helps limit staff exposure to surgical fluids by containing waste throughout collection and transport.
- **Mobile Cart Design:** Smooth-rolling platform allows easy movement between operating rooms and docking areas.
- **Streamlined Disposal Process:** Designed to simplify fluid disposal and reduce manual handling steps.
- **High-Volume Performance:** Built to support demanding surgical environments.



DIMENSIONS

Height: 55 in (140 cm)

Width: 24 in (61 cm)

Depth: 24 in (61 cm)

Weight (empty): 195 lb (88 kg)

Weight (full): 268 lb (122 kg)

GENERAL

Max Vacuum Level: 0-700 mmHg (If connected to hospitable suction, vacuum can be higher than 700 mmHg)

Fluid Capacity: 33,000 mL

Cleaning: Dornoch enzyme solution and Clorox bleach are used to clean the collection cylinders.