

SKYTRON 3008 CLARITY



ABOUT

The **Skytron 3008 Clarity** vascular imaging table is built to support high-quality imaging while providing the flexibility clinicians need during complex vascular and interventional procedures. Designed around an extended carbon fiber tabletop and smooth floating top control, the 3008 Clarity allows precise patient positioning without interrupting workflow. Its fully radiolucent surface helps optimize image quality while supporting lower imaging doses, making it well suited for procedures that rely heavily on fluoroscopy. With intuitive controls and a wide range of motion, this table adapts easily to changing clinical needs across multiple specialties.



FEATURES

- **Extended Radiolucent Carbon Fiber Top:** Features an 83-inch tabletop designed to support clear imaging with minimal interference.
- **Floating Top Functionality:** A single-button float mode allows smooth, multidirectional tabletop movement, even with higher patient loads.
- **Large Imaging Window:** Provides an unobstructed longitudinal imaging area, allowing greater access for vascular and interventional procedures.
- **Intuitive Hand Control:** Ergonomically designed hand control enables easy articulation of the table and floating top.
- **Versatile Clinical Use:** The wide range of motion and imaging access make the 3008 Clarity suitable for vascular, interventional, and imaging-intensive surgical applications.

DIMENSIONS

Table Top Length: 82.5 in (210 cm)

Table Top With: 19.63 in (50 cm)

Width w/ Side Rails: 22 in (56 cm)

ARTICULATIONS

Height Max: 41.75 in (106.5 cm)

Height Low: 27.75 in (71 cm)

Lateral Tilt (left/right): 20°

Trendelenburg: 20°

Reverse Trendelenburg: 20°

Slide X, Y: 8 x 28 in (20.3 x 71.1 cm)

WEIGHT CAPACITIES

Weight Cap Lift & Top Slide: 500 lb (226.8 kg)

Weight Cap Art: 300 lb (136 kg)