

# INTENSA 420



## ABOUT

The **Intensa 420** Procedure and Exam Table is designed to support patients during a wide range of medical examinations while providing clinicians with convenient access and positioning options. Built with an all-metal base, the table offers long-term durability and stability for everyday use in exam rooms and clinical settings. It features multiple storage drawers, including pass-through access, allowing staff to keep essential supplies within reach during procedures. The adjustable components, including stirrups and pelvic tilt, help position patients comfortably while improving access for examinations. Available configurations such as the 420-DW and 420-EO offer additional options to meet specific clinical preferences and workflow needs.



## FEATURES

- **Five Drawer Storage System:** Includes three pass-through drawers and two storage drawers.
- **Adjustable Stirrups:** Five-position continuous-lock stirrups support patient positioning during exams and procedures.
- **Pull-Out Leg Rest with Tray:** Integrated leg rest extends for patient comfort and includes a treatment tray.
- **Durable Metal Construction:** All-metal base provides strength and long-term reliability.
- **Integrated Clinical Features:** Includes pelvic tilt, paper roll holder, and pull-out step.

## DIMENSIONS

**Table Height:** 32 in (81 cm)

**Overall Bed Length:** 60 in (152 cm)

**Overall Bed Length with Leg Extension:** 72 in (183 cm)

**Model 420 Weight:** 275 lb (125 kg)

**Model 420-DW Weight:** 280 lb (127 kg)

**Model 420-EO Weight:** 280 lb (127 kg)

## PATIENT CAPACITY

**Load Capacity (Evenly Distributed):** 500 lb (227 kg)

**Legrest Support Capacity:** 50 lb (23 kg)

**Foot Step Support Capacity:** 800 lb (363 kg)

# SPECIFICATIONS



## TABLETOP AND MOVEMENT

**Backrest Angle Range:** 0° to ~70°

**Upholstered Top Size (W x L):** 27 in x 60 in (69 cm x 152 cm)

**Leg Support Pad Size (W x L):** 18 in x 19.5 in (46 cm x 49.5 cm)

**COM (Customer's Own Material) Requirement:** 3 yd (2.74 m)

## STORAGE

**Pass-Through Drawers (3):** 18 in x 20.5 in x 4 in (46 cm x 52 cm x 10 cm)

**Storage Drawers (2):** 17 in x 17.5 in x 4 in (43 cm x 44 cm x 10 cm)

## ELECTRICAL REQUIREMENTS

**Model 420:** N/A

**Model 420-DW / 420-EO:** 115–250 VAC

## POWER CORD

### (420-DW / 420-EO ONLY)

**Length:** 90 in (229 cm)

**Specifications:** 300V / 18 AWG / 3 conductor / SJT jacket with hospital-grade plug

## ADDITIONAL FEATURES

**Optional Drawer Warmer Temperature:** Approximately 97°F – 99°F

**Paper Roll Sizes (W x Diameter):**

- 18 in x 3 in (46 cm x 8 cm)
- 19 in x 3 in (48 cm x 8 cm)
- 21 in x 3 in (53 cm x 8 cm)

## ENVIRONMENTAL CONDITIONS

**Operating Temperature:** 65°F to 85°F (18°C to 29°C)

**Operating Humidity:** 10% to 90%

**Storage Temperature:** -20°F to 150°F (-29°C to 66°C)

**Storage Humidity:** 10% to 90%