

# ARCOMA MEDSTONE-ELITE 0047-TM5



## ABOUT

---

The **Arcoma Medstone Elite 0047-TM5** pain management table is designed to support patients during pain management, imaging, and fluoroscopy-guided procedures. Built with a radiolucent tabletop and multiple powered movements, the table allows clinicians to precisely position patients while maintaining imaging access throughout the procedure. The TM5 configuration offers five motorized movements, providing flexibility for a variety of interventional pain management and diagnostic applications. Its mobile design and adjustable positioning capabilities help improve workflow while supporting patient comfort and procedural efficiency. With a durable construction and high patient weight capacity, the 0047-TM5 is well suited for pain clinics, surgical centers, and imaging environments.

## FEATURES

---

**Five Movements:** Provides powered height adjustment, Trendelenburg positioning, lateral tilt, longitudinal travel, and lateral tabletop travel.

**Radiolucent Imaging Surface:** Designed to support C-arm and fluoroscopic imaging procedures with unobstructed access to the treatment area.

**Extended Table Travel:** Features both longitudinal and lateral tabletop movement to help clinicians position patients without repeatedly repositioning.

**High Weight Capacity:** Supports patients up to 500 lbs (227 kg) when weight is evenly distributed across the tabletop.

**Mobile Procedure Platform:** Locking caster system allows the table to be easily positioned while providing stability.



## TABLE TOP SIZE

---

**Length:** 82 in (208.3 cm)

**Width:** 22 in (55.9 cm)

## TABLE MOVEMENTS

---

**Lateral Rotation:**  $\pm 15^\circ$

**Vertical Movement:** 11.8 in (30 cm)

**Trendelenburg:**  $\pm 15^\circ$

**Longitudinal Travel:** 15.7 in (40 cm)

**Lateral Travel:** 7.8 in (20 cm)

## WEIGHT

---

**Patient Weight Limit:** 500 lb (227 kg)

**Table Weight:** 650 lb (295 kg)