

# SONY UP-D898DC



## ABOUT

The **Sony UP-D898DC** medical printer is a compact black and white digital thermal printer designed to produce high-quality hard copy images from compatible medical imaging systems. Commonly used with C-arms, ultrasound systems, and cardiac catheterization equipment, the printer delivers fast, detailed prints for patient records, referrals, and clinical documentation. The UP-D898DC features a space-saving design that integrates easily into medical carts and equipment workstations while offering an intuitive LCD display and joystick controls for simple operation. Advanced thermal printing technology produces high-resolution grayscale images without the need for ink or toner, helping reduce routine maintenance. Designed with a DC power input, the UP-D898DC is well suited for mobile imaging applications where DC-powered systems are preferred.



## DIMENSIONS

**Height:** 3.5 in (8.8 cm)

**Width:** 6.2 in (15.4 cm)

**Depth:** 6.5 in (16.5 cm)

**Weight:** 3.9 lb (1.8 kg)

## FEATURES

- **High-Resolution Thermal Printing:** Produces sharp black and white medical images with up to 325 dpi resolution and 256 levels of grayscale.
- **Fast Print Performance:** Generates A6-size prints in approximately two seconds to support efficient clinical workflow.
- **Compact DC-Powered Design:** Space-saving footprint and DC input make the printer ideal for mobile carts and integrated imaging systems.
- **Intuitive User Interface:** Large LCD display and joystick navigation provide quick access to printer settings and status information.
- **USB Image Storage:** Supports saving printed images to a connected USB flash drive for convenient record keeping and image management.

## PRINTER

**Printing Method:** Direct thermal printing

**Resolution:** 325 dpi

**Gradations:** 8 bits (256 levels) of processing

**Picture Elements:** Digital: 4,096 x 1,280 dots

**Picture Area:** 12.625 x 4 in (320 x 100 mm)

**Paper Size:** Width: 4.375 in (110 mm)

### Printing Time

- High Speed Mode: 1.9 seconds / image (960 x 1,280 dots)
- Normal Speed Mode: 3.3 seconds / images (960 x 1,280 dots)

**Picture Memory:** Digital: 4,096 x 1,280 x 8 (bit)

**Interface:** Hi-Speed USB (USB 2.0) (x1)

**Printer Driver Software:** Microsoft Windows XP 32/64 bit, Vista 32/64 bit, 7 32/64 bit, 8 32/64 bit