

BAXTER COMPELLA BARIATRIC BED



ABOUT

The **Baxter Compella Bariatric Bed** is a specialty bed designed for bariatric patient care with a weight capacity of 1,000 lbs. The bed features powered width expansion from 40 to 50 inches and the FlexAfoot bed length adjustment system. It includes integrated caregiver and patient controls, an in-bed scale, and battery backup for continuous operation during power interruptions. The bed supports low positioning with seat deflation to facilitate patient mobility and safety.

FEATURES

- **Powered Width Expansion:** Expands from 40 to 50 inches to accommodate varying patient sizes.
- **FlexAfoot Bed Length Adjustment:** Adjusts bed length for patient comfort and fit.
- **Single-Mode Bed Exit:** Provides bed exit monitoring for patient safety.
- **Digital HOB Angle Indicators and Alert:** Displays head-of-bed angle with alert functionality.
- **Integrated Caregiver and Patient Controls:** Allows operation from both sides of the bed.
- **In-Bed Scale:** Provides patient weight measurement with $\pm 1\%$ accuracy.
- **Battery Backup:** Ensures bed functionality during power loss.



SPECIFICATIONS



DIMENSIONS

Minimum Overall Width: 43 in (109 cm)
Minimum Overall Length: 91 in (231 cm)
Surface Width: 40 to 50 in (102 to 127 cm)
Surface Length (inflated): 80 to 88 in (203 to 224 cm)
Surface Height – Inflated (foam): 8 in (20 cm)
Surface Height – Inflated (low air loss): 8.5 in (21.5 cm)
Sleep Deck to Floor – Low Position: 18.5 in (47 cm)
Sleep Deck to Floor – High Position: 26.8 in (68 cm)
Sleep Deck to Floor – Transport Position: 17 in (43 cm)

BED ANGLES

Head Section: 0 to 50°
Knee Section: 0 to 30°
Max Tilt/Reverse Tilt: 0 to 9°

WEIGHT CAPACITY

Maximum Patient Weight Capacity: 1,000 lbs (453.6 kg)
Scale Weight Capacity: 1,000 lbs (453.6 kg)
Scale Accuracy: $\pm 1\%$ of patient weight (with sleep deck flat)
Weigh in Any Position: Yes

ELECTRICAL INFORMATION

Voltage: 120/230 V
Current: 12/6 A
Frequency: 50/60 Hz

CASTER AND BRAKING SYSTEM

Caster Size: 6 in (15 cm)
Four-Wheel Braking: Yes
Brake Not Set Alert: Yes
Brakes Lock Both Rolling and Swiveling: Yes

POWER AND BATTERY

Battery Backup: Yes

ACCESSORY OPTIONS

IntelliDrive XL Powered Transport: Standard on rental
3-Mode Bed Exit: Standard on rental
Patient-Helper Trapeze: Optional
Patient Pendant: Optional
Smart Bed Ready: Standard on rental
Nurse Call: Capital purchase only
Low Air Loss Surface with Turn Assist: Optional
Low Air Loss Surface with CLRT: Optional