

ENTHERMICS EC250L



ABOUT

The **Entthermics EC250L** is a point-of-care fluid warming cabinet designed for IV and irrigation fluids in clinical settings. The unit features a stainless steel construction with a 2.5 cubic foot capacity, holding up to sixteen 1-liter bags or twelve 1-liter bottles. WarmRight technology ensures fluids are warmed to recommended temperatures for healthcare compliance. The countertop cabinet is stackable and available with optional casters for mobility.

FEATURES

- **WarmRight Technology:** Ensures each bag is safely warmed to the right temperature for IV and irrigation fluids to meet healthcare guidelines.
- **WarmSafe Zone Heating:** Provides evenly balanced heat throughout the chamber for safe and consistent fluid warming.
- **Stainless Steel Construction:** Durable, rugged cabinet designed for long-term reliability in clinical environments.
- **Capacity:** Holds sixteen 1-liter bags or twelve 1-liter bottles in a 2.5 cubic foot compartment.
- **Countertop Desizfiguration:** for flexible placement in patient or procedure rooms.
- **Glide-Out Baskets:** Provides organized storage and easy access to fluids.
- **Optional Casters:** Adds mobility for transport between clinical areas.



SPECIFICATIONS



DIMENSIONS

Overall Height: 22 in (55.9 cm)
Overall Width: 18.5 in (47 cm)
Overall Depth: 23.5 in (59.7 cm)
Compartment Size: 2.5 cubic ft
Compartment Width: 15.5 in (39.4 cm)
Compartment Depth: 18.75 in (47.6 cm)
Compartment Height: 15.25 in (38.7 cm)

TEMPERATURE

Temperature Range (Blanket): 90°F to 160°F (32.2°C to 71.1°C)
Fluid Warming Temperature: Up to 104°F (40°C)

CAPACITY

Total Capacity: 2.5 cubic ft
Bag Capacity: 16 x 1-liter bags
Bottle Capacity: 12 x 1-liter bottles

CONSTRUCTION

Material: Stainless Steel
Shelves/Baskets: 1 basket standard
Optional Glide-Out Baskets: Available
Cooling System: No fans or lint traps

MOUNTING AND MOBILITY

Configuration: Countertop model
Stackable: Yes
Optional Casters: Available

COMPLIANCE

510K Cleared: Yes