

PEDIGO P-2030-S



ABOUT

The **Pedigo P-2030-S** blanket warmer is designed to safely warm and store blankets for hospitals, surgery centers, and other healthcare facilities. Its single-chamber design uses a low-heat-density electrothermal heating system that provides consistent warmth without the need for a circulating fan. Constructed with a durable stainless steel exterior, the cabinet is built for everyday clinical use while offering a simple digital control panel for easy temperature adjustment and monitoring. The spacious interior of the P-2030-S accommodates up to 20 blankets.

FEATURES

- **Consistent Warming:** Low-heat-density electrothermal heating system provides even heat distribution without a circulation fan.
- **Digital Control:** User-friendly control panel features a digital display with adjustable temperatures up to 180°F (83°C).
- **Durable Construction:** Stainless steel exterior is designed to withstand daily use.
- **Efficient Storage:** Interior capacity holds approximately 15 to 20 blankets.



DIMENSIONS

Height: 40.4 in (102.5 cm)

Width: 22.3 in (56.7 cm)

Depth: 31.6 in (80.3 cm)

Weight: 175 lb (79 kg)

HEATING CONTROLS

Adjustable control with digital display

Temperature Range: Adjustable up to 180°F (83°C)

INTERIOR DIMENSIONS

Interior Capacity: 7.7 cubic feet (0.2 m³)

Blanket Capacity: 15-20 blankets