

PEDIGO P-2040



ABOUT

The **Pedigo P-2040** blanket warmer is designed to safely warm and store blankets for high-volume patient care areas, helping improve patient comfort before, during, and after procedures. Its dual-chamber design provides increased storage capacity while allowing staff to organize blankets for efficient access throughout the day. Using a low-heat-density electrothermal heating system, the cabinet delivers consistent, even warmth without the need for circulating fans, reducing maintenance and noise. Independent digital controls with programmable timer functions allow users to manage temperatures. Built with a durable stainless steel exterior and mounted on heavy-duty casters, the P-2040 is well suited for operating rooms, recovery departments, labor and delivery units, and other demanding clinical environments.



FEATURES

- **Dual-Chamber Storage:** Two warming compartments provide increased capacity for facilities with higher blanket usage.
- **Even Heating:** Low-heat-density electrothermal heating system delivers consistent blanket warming without circulating fans.
- **Programmable Controls:** LED temperature display, adjustable settings, timer function, and lockout features.
- **Safety Features:** Includes an independent over-temperature protection system and safety override to help prevent overheating.
- **Durable Design:** Stainless steel construction, insulated cabinet, heavy-duty locking casters, and epoxy-coated interior components.

DIMENSIONS

Height: 73.2 in (185.8 cm)

Width: 22.1 in (56.1 cm)

Depth: 32.3 in (82.7 cm)

INTERIOR DIMENSIONS

Upper Compartment: 7.7 cubic ft; 26.75 x 19 x 26.5 in (HxWxD)

Lower Compartment: 7.7 cubic ft; 26.75 x 19 x 26.25 in (HxWxD)

Total Capacity: 60 one-liter bottles or 40-50 blankets

HEATING CONTROLS

Two compartments, each with individual controls

Blanket Temp Range: 90°F – 160°F

Irrigation Temp Range: 90°F – 150°F