

STERIS 4095



ABOUT

The **Steris 4095** general surgical table is designed to provide flexible patient positioning, imaging access, and reliable support across a wide range of surgical specialties. Built with a radiolucent tabletop and a long imaging window, the table is well suited for procedures that require frequent C-arm access, including orthopedic, neurosurgical, and cardiothoracic cases. Its powered tabletop slide and articulation features allow the surgical team to position the patient without repeatedly reconfiguring the table, helping improve efficiency in the operating room. The dual-articulating headrest adds versatility for anesthesia access and airway management, while the modular design supports a variety of surgical setups.

FEATURES

- **Radiolucent Imaging:** Large imaging window and radiolucent tabletop help provide clear C-arm access.
- **Tabletop Slide:** Offers extended longitudinal slide for improved patient positioning and surgical access.
- **Dual-Articulating Headrest:** Supports flexible head and upper body positioning for anesthesia access, airway management, and specialty procedures.
- **Versatile Positioning:** Designed to accommodate a wide range of positions including supine, prone, lateral, Fowler's, and more.
- **Mobile, Modular Design:** Features a maneuverable base and removable leg section.



DIMENSIONS

Height Range: 26 to 45 in (67.5 to 113.6 cm)

Width: 20 in (50.8 cm)

Width with Side Rails: 22 in (54.7 cm)

Length: 81 in (205.7 cm)

Shipping Weight: 560 lb (254 kg)

PATIENT WEIGHT CAPACITY

1,100 lb (499 kg): Patient support capacity, including raise/lower, centered on the column

1,000 lb (454 kg): Full table articulation, centered on the column

600 lb (272 kg): Full table articulation including slide in Normal Orientation, and slide only toward the feet in Reverse Orientation

575 lb (261 kg): Full table articulation including slide toward the head in Reverse Orientation

SPECIFICATIONS



TABLE MOTION & POSITIONING

Table Slide Range: 9 in (22.7 cm) toward the head, 9 in (22.7 cm) toward the foot

Trendelenburg / Reverse Trendelenburg: 30° / 30°

Lateral Tilt (Left / Right): 20° / 20°

Flex / Reflex: 140° / 100°

ADDITIONAL FEATURES

Manual Override

Perineal Cut-Out

Kidney Elevator

SECTION ARTICULATION

Head Section: +90° / -90°

Back Section: +80° / -40°

Leg Section: 0° / -105° (removable)