

STRYKER 1000



ABOUT

The **Stryker 1000** is a manual emergency care stretcher that uses a hydraulic pump base and a pneumatic-assisted Fowler section for positioning. Height is adjusted with foot pedals connected to a hydraulic jack at each end of the base, moving the stretcher through a range of roughly 22" to 35". The Fowler backrest adjusts from flat to 90 degrees, the knee gatch section raises up to 35 degrees, and the base can be tilted into Trendelenburg or reverse Trendelenburg between +12 and -18 degrees. It carries a rated weight capacity of 500 pounds on a 78" x 28" bed surface.



FEATURES

- **Hydraulic Height Adjustment:** Foot pedal-operated hydraulic jacks raise and lower the stretcher between approximately 22" and 35".
- **Pneumatic Fowler Section:** Backrest adjusts from flat to 90 degrees using a gas-assisted Fowler mechanism.
- **Knee Gatch and Trendelenburg:** Knee gatch adjusts up to 35 degrees, with base tilt from +12 to -18 degrees for Trendelenburg and reverse Trendelenburg.
- **500-Pound Weight Capacity:** Rated for a maximum patient weight of 500 pounds on a 78" x 28" bed surface.

OVERALL DIMENSIONS

Overall Length: 78 in (198.1 cm)

Overall Width: 28 in (71.1 cm)

Height Range: 22 in (55.9 cm) - 35 in (88.9 cm)

WEIGHT CAPACITY

Maximum Weight Capacity: 500 lb (226.8 kg)

POSITIONING

Fowler Angle: 0° - 90°

Knee Gatch Angle: 0° - 35°

Trendelenburg / Reverse: +12° to -18°