

STRYKER SURGISTOOL

ABOUT

The **Stryker Surgistool** is a surgeon's stool designed for use during surgical procedures, with adjustable seating to accommodate different positioning preferences. It is available with a teardrop or wedge seat, and the contoured backrest adjusts to provide additional support. A foot-controlled hydraulic base raises and lowers the seat, allowing height adjustments without using your hands. The three-spindle base provides easy access to foot pedals while maintaining a stable, mobile design. It also features a front-end locking system to help keep the stool in place when needed.

FEATURES

Seat Options: Available with a teardrop or wedge seat to accommodate different positioning preferences.

Adjustable Backrest: Contoured backrest adjusts to provide support and help the user find a comfortable sitting position.

Dual-Density Foam Cushion: Foam cushioning combines firm and soft layers for support and comfort.

Foot-Controlled Height: Hydraulic adjustment provides a height range of 21 to 28.5 inches, or 23 to 30.5 inches with the optional riser kit.

300 lb Capacity: Supports users weighing up to 300 lbs (136 kg).



SPECIFICATIONS



DIMENSIONS

Height Range: 21 to 28.5 in (53 to 72 cm)

Height Range (with Riser Kit): 23 to 30.5 in (58 to 77 cm)

Weight Capacity: 300 lb (136 kg)

SEATING

Seat Style: Teardrop or wedge

Seat Cushion: Dual-density foam

Backrest: Adjustable, contoured

Armrests: Multi-position padded armrests (Optional)

BASE & MOBILITY

Base: Three-spindle base for stability and easy pedal access

Height Adjustment: Foot-controlled hydraulics

Locking: Front-end locking system

Hood: Thermoformed plastic hood

OPTIONS

Riser Kit: 2 in (5 cm), raises the seat height range (Optional)

Armrests: Multi-position padded armrests (Optional)