

STERIS MEDIVATORS SCOPE BUDDY PLUS



ABOUT

The **Steris Medivators Scope Buddy Plus** is an endoscope flushing aid. It pumps fluid through endoscope channels the same way a technician would with a syringe, but automatically, through a color touchscreen with manual or automatic mode. The system doses detergent on its own and tracks each flushing cycle, so staff can pull up records by user, endoscope, and patient ID with the optional barcode scanner. It works with most current flexible endoscopes, including those with elevator wire channels, and the 24-hour disposable tubing cuts down on the contamination risk that comes with reusable lines. Like the original Scope Buddy, it's not a disinfectant, so scopes still need high-level disinfection afterward.



FEATURES

Automated Detergent Dosing: Doses chemistry automatically during flushing instead of requiring manual measuring.

Manual or Automatic Mode: Both modes walk staff through a consistent, step-by-step flushing process.

Track and Trace Logging: Records user, endoscope, and patient IDs for each cycle, exportable for reporting.

24-Hour Disposable Tubing: Reduces contamination risk compared to reusable tubing lines.

Broad Endoscope Compatibility: Works with most current flexible endoscopes, including elevator wire channel models.



SPECIFICATIONS



DIMENSIONS

Height: 10.3 in (26.2 cm)

Width: 7.5 in (19 cm)

Depth: 11.5 in (29.2 cm)

Weight: 10 lb (4.5 kg)

OPERATION

Maximum Pump Pressure: 23 psi (1.6 bar)

Dosing Accuracy: ± 2 mL (for 0 – 20 mL doses); $\pm 10\%$ (for doses over 20 mL)

Electrical Rating: 100 – 240 V, 50/60 Hz

Power Supply: Cord connection

Overvoltage Category: II

Pollution Degree: 2

ENVIRONMENTAL SPECIFICATIONS

Operating Temperature: 59°F to 95°F (15°C to 35°C)

Operating Humidity: 10 – 75% (non-condensing)

Storage Temperature: -20.2°F to 134.6°F (-29°C to 57°C)

Storage Humidity: 5 – 90% (non-condensing)

Altitude: Under 6,526 ft (2,000 m)